

Class:3rd stage

Subject: Thermodynamics

م.م زهراء عبد الاله هادي م.م ايلاف جاسم موهان م.م حسن علي شمخي مهندس احمد عباس محمد



Ministry of Higher Education and Scientific Research

Al-Mustaqbal University College Chemical Engineering and Petroleum Industries (Thermodynamic Lab3)

Experiment No. 1
(Air speed measurement)

اعداد:
م.م زهراء عبد الاله هادي
م.م ايلاف جاسم موهان
م.م حسن علي شمخي
مهندس احمد عباس محمد



Class :3rd stage
Subject: Thermodynamics

م.م زهراء عبد الآله هادي م.م ایلاف جاسم موهان م.م حسن علي شمخي مهندس احمد عباس محمد



Number of Experiment: Experiment number one.

Name of Experiment: Air speed measurement.

Purpose: Calculate air speed.

Equipment: Used equipment in this experiment are:

- 1. An emometer which is used to measure air speed (V_{air}) as shown in Figure below.
- 2. Data Logar which has two thermocouple and used to measure inlet and outlet temperature.
- 3. Blower and control unit.



Fig: Anemometer.



Class:3rd stage

Subject: Thermodynamics

م م زهراء عبد الآله هادي م م ايلاف جاسم موهان م م حسن علي شمخي مهندس احمد عباس محمد



Theory:

Air is all around us, but we cannot see it. Gravity from the Earth pulls air down - this is called air pressure. We don't feel this pressure because our bodies push an equal amount of pressure outward.

Everything is made from molecules, the air also contains molecules. The temperature outside on any given day is determined by the speed of air molecules.

When the airspeed increases the temperature increases too due to the increase of air molecules movement and vise versa. Generally, the larger the temperature difference, the larger the resulting airspeed will be.

The airspeed is important information for accurate navigation of an aircraft. To maintain a desired ground track whilst flying in a moving air mass, the pilot of an aircraft must use knowledge of wind speed, wind direction, and true air speed to determine the required heading. So it is important to measure the airspeed and this done by an instrument called an anemometer which measures the airspeed by monitoring the amount of heat removed from a surface using one or more simple temperature sensors.



Class:3rd stage

Subject: Thermodynamics

م.م زهراء عبد الآله هادي م.م ایلاف جاسم موهان م.م حسن علي شمخي مهندس احمد عباس محمد



Procedure:

- 1. Fix an emometer and the two thermocouple (one for the inlet temperature (T_{in}) and the other for the outlet temperature (T_{out})).
- 2. Record outlet temperature which is constant in the experiment.
- 3. Take the first reading of speed and inlet temperature.
- 4. After 1-2 minutes, change the speed and record T_{in} and V again and repeat this step until reach the required number of readings.
- 5. Make a table with results as below:

No.	T _{out} (°C)	T _{in} (°C)	V _{air} (m/s)
1			
2			
3			

Requirement:

- 1. Draw curve between T_{in} and V_{air} at constant T_{out} .
- 2. Define anemometer.
- 3. Calculate mach number (M_a) by using equation below:

$$M_{a=}\frac{\textit{V}}{661.9*\sqrt{\frac{\textit{T(in)}}{\textit{T(out)}}}}$$