

Practical Pharmacognosy

3rd. Stage

1st semester

Dr. Zahraa Shubber

Lab.3




Tannins

- Tannins comprise a large group of complex substances that are widely distributed in the plant kingdom.
- Chemically tannins are **complex substances**; they usually occur as mixtures of **poly hydroxyl phenols** that are difficult to separate because they **do not crystallize**.
- Tannins are divided according to the identity of the phenolic nuclei involved, and on the way they are joined into two classes:




1 .Hydrolysable tannins:

- This class consists of **gallic acid** and related polyhydroxy compounds (**hexahydroxydiphenic acid**) and their derivatives esterified with glucose.
- They are termed **hydrolysable tannins** due to **ease of esters to hydrolyze** to phenolic acids and sugar.
- They were formerly known as **pyrogallol** tannins.




Gallic acid




Hexahydroxydiphenic acid

2. Nonhydrolysable tannins or condensed tannins:

- This class contains **only phenolic nuclei** but frequently linked to **carbohydrates or proteins**.
- When treated with hydrolytic agents, these tannins tend to polymerize, yield insoluble usually **red-colored** products known as **phlobaphenes**.
- The name (condensed tannins) is due to the fact that on the treatment with hot acid some of C-C bonds are broken yielding **anthocyanidin monomers**. These tannins are sometimes called **catechol tannins**.



Catechin



Anthocyanins

(Are glucosides of anthocyanidins)

General properties of Tannins:

- Tannins are non-crystallizable compounds that, with *water* form *colloidal solution* possessing acid reaction and sharp "*puckering*" taste.
- They cause precipitation of solution of *gelatin* as well as *alkaloids*.
- They form dark *blue, greenish* black soluble compounds with *ferric salts*.
- They produce deep *red* color with *potassium ferricyanide* and *ammonia*.
- They are precipitated by salts of *copper, lead* and *tin* by strong aqueous *potassium dichromate* or *1% chromic acid* solution. In alkaline solutions; many of their derivatives readily absorb oxygen.
- Tannins precipitate *proteins* from solution and can combine with proteins, rendering them resistant to proteolytic enzymes. When applied to living tissue this action is known as an "*astringent*" action and form the basis for therapeutic application of tannins.

Uses of tannins:

1. **Astringents**, used in the gastrointestinal tract and on the skin abrasion.
2. In the **treatment of burns**, the proteins of the exposed tissue are precipitated and form a mildly antiseptic protective coat under which the regeneration of new tissue may take place.
3. Use in the process of vegetable- tanning which converts animal hides to leather (**leather industry**).
4. **Antidote** treatment of alkaloids poisoning.
5. **Ink industry**.



Pyrogallotannins(Nut gall)



Catechole tannins(Hamamelis leaf)

Isolation & Identification of the Tannins.

Aim: Isolation of the Tannins.

Equipment :

- Medium size beaker.
- Conical flasks.
- Heater.
- Filter paper.
- Funnels.

Reagents:

Distill water.

Procedure

Method of extraction: Decoction.

Plant used: Cinnamon.

family: Lauraceae.

Part used: Dry bark.

Add *2gm* of cinnamon, coarsely powdered to a beaker.



Add *30ml* of distill water.



Boil for *5-10 mins*.



Filter.



The filtrate will resemble the *Tannin* extract.

The Chemical Tests:

(Ferric Chloride Test)

Aim: Identity test for Catechol Tannins.

Equipment & Reagents:

- Test tube.
- Pipette.
- Ferric chloride.
- Distill water.

Procedure:

- Add 2ml of the extract to the test tube.
- Add 2ml of distill water.
- Add 2 drops of ferric chloride.

Results:

Green-Black color will be produced which indicates condensed tannins.



Discussion

- The ferric chloride test is used to determine the presence of the **phenolic nuclei** present in the tannin molecule.
- Phenols form a complex with ferric ions. This complex has an intense dark-green color.
- $6 \text{ PhOH} + \text{Fe}^{3+} \rightarrow [\text{Fe}(\text{OPh})_6]^{3-}$



Thank You