



**Ministry of Higher  
Education and  
Scientific Research**

**AL Mustaqbal university college  
Chemical Engineering Department**

**Analytical chemistry laboratory**

**Experimental No.1  
Preparation of standard solution of  
 $\text{Na}_2\text{CO}_3$**

**Prepared by**

**Asst. lect. Ban Ali Hassan**

**2021-2022**

- **Purpose of This Test:**

**Determination Normality of Solution.**

- **Introduction and theory :**

- **Primary Standard solution:**

The primary standard is a highly purified compound that serves as a reference material in volumetric and mass titrimetric methods. The accuracy of a method is critically dependent on the properties of this compound. Important requirements for primary standard are the following: High purity (e.g.99.9%by weight), atmospheric stability, absent of hydrate water, modest cost, reasonable solubility in the titration medium and reasonable large molar mass. One of such compounds is sodium carbonate ( $\text{Na}_2\text{CO}_3$ ).

- **Secondary Standard solution:**

Some compounds consider as secondary standard solution due missing one of the important requirements for primary standard. For that can prepare as an approximate normal solution and then calibrate with a known standard chemical solution by volumetric analysis. Such as: hydrochloric acid HCL, sodium hydroxide NaOH.

## Procedure :

- To prepare 250ml of 0.1N Sodium carbonate (Na<sub>2</sub>CO<sub>3</sub>).

- 1- Dry Appropriate quantity of sodium carbonate in an oven set at 110°C for 1hr. Store the dried material in desiccators.
- 2- Calculate the weight of (Na<sub>2</sub>CO<sub>3</sub>) from the law.
- 3- Weigh (Na<sub>2</sub>CO<sub>3</sub>) using a microbalance (1.0 or 0.1) mg accuracy.
- 4- Transfer this quantity into a 250 ml beaker and solve it in about 50ml distilled water (D.W). Transfer the solution into a 250ml volumetric flask. Complete the solution to the mark by distilled water with the same beaker. Stopper the volumetric flask and shake well.

## Calculation :

$$N = \frac{\text{wt}}{\text{eq.wt}} \times \frac{1000}{V \text{ in ml}} \text{ eq /L}$$

$$\text{eq.wt} = \frac{\text{M.wt}}{\text{No. of eq}}$$

W<sub>t</sub> = Weight of Na<sub>2</sub>CO<sub>3</sub>(gm)

V = volume of volumetric flask (ml)

N = Normality of solution.

## Discussion:

- Determine the volume (ml) required to dilute a solution containing 48.480g of K<sub>2</sub>SO<sub>4</sub> in order to prepare a 1 N of the solution.
- What are the specifications of primary standard solution?

• **Preparation of approximately (0.1N) HCl from 37% and 1.2 Sp.g**

1- Calculate the Normality of the concentrated HCl from the law:

$$N_1 = \frac{\rho \times \% \times 10}{\text{eq.wt.}}$$

$$\text{Sp. Gr.} = \frac{\rho_{\text{HCl}}}{\rho_{\text{water}}}$$

$$\text{Eq. wt.} = \frac{\text{Mwt. of HCl}}{\text{no of H}^+}$$

2- To prepare (500ml) of 0.1N HCL:

$$N_1 V_1 = N_2 V_2$$
$$N_1 \times ? = 0.1 \times 500$$

3- Transfer some of distilled water to 500ml volumetric flask and add the calculated volume of conc. HCl. Mix the content thoroughly. Complete to the mark with distilled water and make a label.

- **Discussion:**

- What are the specifications of secondary standard solution?
- Determine the volume (ml) required to dilute a solution of HCl from conc. 37% and Sp.g. 1.2 to 1000ml if the Normality is 0.25 eq/l