

ANALYTICAL CHEMISTRY LAB

Determination of Chloride by Mohr method


The first stage-First semester

M.s.c. Zaynab Faleh




Precipitation Titration

Mohr's method used to determine the chloride ion concentration of a solution by titration with silver nitrate ($AgNO_3$). A soluble chromate salt (K_2CrO_4) is added as the indicator. This produces a yellow color solution as the silver nitrate solution is slowly added, a precipitate of silver chloride formed.



This method uses chromate as an indicator, chromate forms a precipitate with Ag ions, but this precipitate has a greater solubility than that of ($AgCl$) therefore ($AgCl$) is formed first and after all (Cl) ions are consumed. The end point of the titration occurs when all the chloride ions are precipitated, then additional silver ions react with the chromate ions of indicator, potassium chromate, to form a red brown precipitate of silver chromate.



Purpose of this experiment

- ▶ Determination the concentration of silver in precipitation method

- ▶ **Chemicals and tools**

NaCl, AgNO₃, K₂CrO₇ as indicator, distilled water, burette, pipette, stand, clamp, brush, conical flask, spatula, funnel, volumetric flask, washing bottle, beaker, dropper, balance, watch glass.

Experimental part

- ▶ 1-Fill the burette with (0.1) N of AgNO_3 .
- ▶ 2- Take (5) ml of NaCl sample by pipette and put it in a conical flask.
- ▶ 3- Add (5) drops of K_2CrO_7 indicator to the to the conical and mix well
- ▶ 4- Titrate with AgNO_3 until the appearance of red-brown precipitate.
- ▶ 5- Repeat the titration 3 times and take the average

Calculation

- ▶ Calculate the concentration of Chloride ion in normality by using the law:

$$(N \times V)_{\text{AgNO}_3} = (N \times V)_{\text{NaCl}}$$

(Average = $V_1 + V_2 + V_3 / 3$)

